

H. pylori RUT TEST KIT

(Rapid Urease Test Kit)

A rapid, reliable and easy-to-use test for the qualitative detection of *Helicobacter pylori* in gastric biopsy specimens.



The Gastro Hub® H. pylori RUT Test Kit is based on the rapid urease reaction, providing clear color-change results to assist clinicians in the diagnosis of H. pylori infection and support timely treatment decisions.

KEY FEATURES

- Rapid urease-based detection of H. pylori
- Easy-to-read color change for result interpretation
- High sensitivity and specificity
- Simple, ready-to-use test format
- Suitable for gastric biopsy specimens collected during endoscopy
- Stable and reliable performance
- Supports early diagnosis and treatment planning
- Manufactured under strict quality standards




RESULT INTERPRETATION

POSITIVE RESULT (+)



The test medium changes from yellow/orange to pink or red, indicating the presence of urease-producing *Helicobacter pylori* in the biopsy specimen.

NEGATIVE RESULT (-)



The test medium remains yellow/orange with no significant color change, indicating that *Helicobacter pylori* was not detected under the test conditions.



Note:

Test results should always be interpreted in conjunction with the patient's clinical findings, endoscopic observations, and other laboratory investigations.

HOW TO USE

- Collect biopsy sample during endoscopy.
- Place the biopsy specimen on the test medium.
- Close the lid and incubate at room temperature.
- Observe color change within 1 – 24 hours.
- Interpret the results based on color change.




14. Hezell S- L., Borody T. J., Gall A. and Lee A. :Campylobacter Pyloridis Gastritis I: detection of urease as a marker of bacterial colonisation and gastritis, Amer J. gastroenterol. 1987, 82, 192-196.
15. Rauws Eric A. J., Weies Langenberg, Henrik J, Houthof H- C. Zanen and GNJ Tytgat: Campylobacter Pyloridis associated active chronic gastritis; a prospective study of its prevalence and the effect of antibacterial and antiulcer treatment: Gastroenterology 1988, 94, 33-40,
16. Deltenre M., Glupczynski Y., Nyst J. F., de Koster E., Burette A., Assessment of four modalities of urease test in the diagnosis and follow-up of Campylobacter Pylori associated gastritis. Gastroenterology : 1989, 96 (5 pt 2)A117.
17. Deltenre M., Glupczynski Y., C de Prez, Nyst J. F., Burette A., Labbe M., Jonas C., de Koster E. The reliability of urease tests, histology and culture in the diagnosis of Campylobacter Pylori infection; ScandJ-Gastroenterol 1989;24(Suppl 160): 19-24.
18. Goodwin C. S., Armstrong J-A., Chilvers T.; transfer of Campylobacter Pylori and Campylobacter Mustlae to Helicobacter Genov as Helicobacter comb nov and Helicobacter Mustlae respectively; Int. J. Syst. Bacterol., 1989, 39, 397-405.
19. Calam J., Helicobacter Pylori, Medicine Digest. 5 May 1990, 5.
20. Tytgat G. N. J., Axon A. T. R., Dixon M. F., Graham D. Y., Lee A, and Marshall B. J.: Helicobacter Pylori : Causal agent in peptic ulcer disease?: World congress of gastroenterology Sydney, Australia. 1990, Working Party report page 36-45.
21. Kang J. Y. Wee. A. : Helicobacter Pylori and gastric acid output in peptic ulcer disease; Digestic diseases and sciences Vol 36 No. 1 January 1991 5-9.
22. Logan R.P.H., Grammett P.A., Misiewicz J.J., Karim Q.M., Walker M. M& Barom J. H.: One week eradication regimen for Helicobacter Pylori : The Lancet, Vol. 338, Nov. 16, 1991, 124-152.
23. Graham D. Y. Lew G. M., Klein P. D. et al; Effect of treatment of Helicobacter Pylori in long term recurrence of gastric and duodenal ulcer; Ann. Int. Med., 1992., 116. 705-708.
24. Borody T, Andrews P, Mancuso N., Jankiewicz E., Brandl S.; Helicobacter reinfection 4 years post eradication: Lancet 1992, 339. 1295.
25. Hentschel E., Brandstatter G., Dragosics B., Hirschl A. M., Nemec H., Schutze K., Tauber M., Wurzer H.; Effect of Ranitidine and Amoxycillin plus Metronidazole on eradication of Helicobacter Pylori and the recurrence of duodenal ulcer: New England Journal of Medicine Feb 4, 1993. Vol. 428 No. 5.
26. Eurogast Study Group. An international association between Helicobacter Pylori infection and gastric cancer. Lancet 1993; 341 : 1359-62.
27. Witherspoon A.C., Doglioni C, Diss TC, Pan K, Moschine A, Boni M, et al. Regression of primary low grade B-cell gastric lymphoma of mucosa-associated lymphoid tissue type after reacidation of Helicobacter Pylori. Lancet 1993; 342: 575-7,
28. Slotte M. Eidt S. Healing gastric (MALT)-lymphomas by eradicating H. Pylori, Lancet 1993; 342: 568.
29. Bayerdorffer E. Eidt S. Lehn N, Thiede D, et al. Double-blind treatment of early gastricmalt-lymphoma patients by H. Pylori eradication (MALT). Gastroenterology 1994; 106: 370.

A DIAGNOSTIC KIT FOR EXTERNAL USE ONLY.


DETECTION OF UREASE ENZYME OF H. PYLORI

PYLO DRY™

• Easy to carry out • Dependable • Economic • Convenient



Biopsy Material Being Inserted into the dry Media & add one Drop of Distilled Water.




Test is Positive for H. Pylori



Store the Kits at Room Temp. & Dry Place

RESPECTIVE & PRINCIPAL



Dr. J.R. Warren

Fitzgerald and Murphy noted an association between urease production in the gastric mucous layer and peptic ulcer disease in 1950. Warren and Marshall carried out evaluation of all patients with dyspeptic symptoms requiring upper G. I. Endoscopy and described a *Campylobacter* like organism in gastric mucosa in 1983. It was first named as *Campylobacter pyloritis*, later changed as *Campylobacter Pylori*. But the observations over the next few years have revealed, this organism has very little resemblance in character with *Campylobacter*. So, in 1989 it was renamed as *Helicobacter*, because of its helical shape.

H. Pylori is associated with duodenal ulcer in 60-80%. *H. Pylori* also plays a significant role in non-ulcer disease- Presently it is considered to be the major cause of the chronic gastritis. Latest studies also show evidence of definite association of *H. Pylori* with lymphoma and carcinoma of the stomach. Remission and exacerbation was the accepted natural history of duodenal ulcer but since the possible role of *H. Pylori* was recognised, previous acid hypothesis alone is not sufficient to explain this fact,



Dr. B.J. Marshall

Bismuth compounds, antibiotics, Nitromidazole and later Omeprazole alone was found effective in eradication of *H. Pylori* but developed high incidence of relapse. However, dual therapy (e.g. Omeprazole & Amoxicillin) and triple therapy with various combinations of antisecretory drugs (e.g. Omeprazole, Ranitidine etc.) and antibiotics (e.g. Amoxicillin, Tetracycline and newer agent Clarithromycin) have very high eradication rates.



After eradication reinfection with *H. Pylori* is rare. The enzyme urease is a high molecular weight protein and very active in comparison to the urease enzyme produced by other bacteria. High concentration of urease is present in gastric mucosa. *H. Pylori* can tolerate a p^H value even below 3.0 due to rapid formation of ammonia which raises the p^H immediately.

Urease reacts with urea in the presence of water molecules and Hydroxide ions producing ammonia.

The presence of *H. Pylori* in the biopsy taken from gastric mucosa of the patients can be done by in vitro urease test which is of course simple and accurate.

COMPOSITION-Ureca, a pH indicator, bacterostatic agent, a stabilising agent.

PROCEDURE OF THE TEST

Biopsy should be taken from the pyloric antrum adjacent to the opening avoiding the areas of erosion and ulceration. Pull up the sticker and introduce the biopsy in the exposed yellow media and add one drop of distilled water.

Place the covering sticker as before.

Urease enzyme of *H. Pylori* reacts with urea of the media and changes the colour from yellow to red or pink elevating pH.

OBSERVATION - Observe the colour change from yellow to red or pink at 10 minutes intervals for four hours. The observation should be continued upto 24 hours. If the colour changes from yellow to red or pink the result is positive.

PRECAUTIONS:

Agents which cause suppression of *H. Pylori* e.g. Antibiotics, Bismuth, Nitromidazole and high doses of antisecretory drugs, should be avoided for 3 weeks prior to this test to prevent false negative results.

Repeat endoscopy and biopsy should be done to check eradication about one month after completion of the regimen.

REFERENCES

1. Krientitz W. Weber das Auftreten von mageninhalt bei Carcinoma Ventruculi. Deutsche Medizinische Wochenschrift 22: 872, 1906.
2. Rosenow, E. C., Sandford, A. H. The bacteriology of ulcer of the stomach and duodenum in man J infect Dis 1915, 17: 210-226.
3. Doenges, J. L. Spirochaetes in the gastric glands of Macacus rhesus and of man without related disease. Arch Pathol 1939, 27: 469-477.
4. Pitzzi Sald O., Murphy P.: Studies on Physiological Chemistry and clinical significance of urease with reference to stomach; Irish J. Med. Sci, 1950, 292: 99-153.
5. Warren J. R. : Unidentified curved bacilli on gastric epithelium in active chronic gastritis, June 1983, The Lancet. Letter, 1, 1273.
6. Marshall B. J.: Unidentified curved bacilli on gastric epithelium in active chronic gastritis, June 1983, The Lancet, Letter, 1, 1273.
7. Langenberg M. L. : Tytgat G. N. J., et al: 1984 Lancet 1; 1348.
8. Morris A., Nicholson G.: Campylobacter Pyloridis; ingestion causing gastritis and acrohydria, Austral Microbiol. 1986, 7, 206.
9. Goodwin C. S. Armstrong J. A., Marshall B.: Campylobacter Pyloridis and peptic ulceration, J. Chile., Pathol; 1986; 39, 533.
10. Waghorn D.: Campylobacter Pyloridis : a new organism to explain an old problem? ; Postgrad- Medical Journal 1987, 63 533.
11. Marshall B. J., Goodwin C. S. 1987 : Revised nomenclature of Campylobacter Pylori. International Journal of systematic Bacteriology, 37:68.
12. Coghlan J. G., Humphrys H., Cooley C., Keane C., Gilligan D. Sweeney E. O'Morain C.: Campylobacter Pylori and recurrence of duodenal ulcer - A 12 month Study. Lancet 1987. Nov. 14, 1109-1111.
13. Marshall B. J., Warren J. R., Francis G. J. et al: Repaid urease test in the management of Campylobacter Pyloridis associated gastritis; Amer. J. Gasroenterol, 1987, 82, 200-210.