

Multi Band Ligator Barrel with 6 loaded Bands.



Used for Endoscopic Ligation of Esophageal Varices at or above Gastro Esophageal Junction.

This device is supplied non sterile, individually packed - intended for single use only.

Technical Information :

Shaili Band Ligator requires 9.5 to 11.2 .mm endoscope diameter.


Shaili Paediatric Band Ligator requires 8.0 to 9.5 mm endoscope diameter.

Band Ligator Set Contains

- Applicator x 1 pc., Barrel loaded with 6 x Bands,
- Barrel loaded with 7 x Bands,
- Atraumatic plastic hook with Tube X1pc.,
- Flushing Needle blunt tip x1 pc.
- Second last band with different colour.

Barrel Contains

- Barrel loaded with 6 x Bands,
- Barrel loaded with 7 x Bands,
- Atraumatic plastic hook with Tube X1pc.



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Reference	Product Name
SBLS-6-F*	Multi Band Ligator Set	SBLS-6-O*	Multi Band Ligator Set
SBLS-7-F*	Multi Band Ligator Set	SBLS-7-O*	Multi Band Ligator Set
SBL-6-F*	Multi Band Ligator Barrel	SBL-6-O*	Multi Band Ligator Barrel
SBL-7-F*	Multi Band Ligator Barrel	SBL-7-O*	Multi Band Ligator Barrel
SBLS-6-P-F*	Multi Band Ligator Set Paediatric	SBLS-6-P-O*	Multi Band Ligator Set Paediatric
SBL-6-P-F*	Multi Band Ligator Barrel Paediatric	SBL-6-P-O*	Multi Band Ligator Barrel Paediatric
SBLS-7-P-F*	Multi Band Ligator Set Paediatric	SBLS-7-P-O*	Multi Band Ligator Set Paediatric
SBL-7-P-F*	Multi Band Ligator Barrel Paediatric	SBL-7-P-O*	Multi Band Ligator Barrel Paediatric

* To Indicate the Length of Thread., F = 175 cm. for Fujinon Endoscope, O = 155 cm. For Olympus / Pentax Endoscope.