

# PEG Kit

Used for percutaneous endoscopic gastrostomy placement  
to provide enteral nutrition to the patient requiring nutritional support

## Technical Information

Size - 20 - 22 - 24 Fr. Pull technique  
Other size can be produced on request.  
Availability : Sterile-for single use only.

## PEG Kit CONTAINS :

- 1 Silicon feeding tube.
- 1 Loop insertion wire
- 1 Cable tie
- 1 Universal adaptor.
- 1 Bolus adaptor.
- 1 Silicon bush for locking.
- Introducer needle - Big needle
- 11 No. scalpel
- 5 ml. Syringe with needle.
- 25 Gauge needle - Small needle.
- 1 Scissor.
- 1 Instruction manual.
- Availability of snare is optional.



## Ordering Information- - Mention product code x quantity required

Reference	Product Name	Size	Remarks
S-PEG-20	PEG Kit	20FR	
S-PEG-20-SN	PEG Kit	20FR	With Snare
S-PEG-22	PEG Kit	22FR	
S-PEG-22-SN	PEG Kit	22FR	With Snare
S-PEG-24	PEG Kit	24FR	
S-PEG-24-SN	PEG Kit	24FR	With Snare