

Naso Jejunal Feeding Tube



Nasojejunal Feeding Tubes are frequently used for long-term enteral nutrition


Technical Information

Feeding tube size diameter - 10, 12, 14, & 16 fr

Working length - 120 cm & 155 cm

Wire guide compatible -0.038 inch.

Sterile, single use only.



Ordering Information : Mention the product code x quantity required.

| Reference | Product Name | Size | Length |
|--------------|---------------------------|------|--------|
| S-NJFT-10120 | Naso Jejunal Feeding Tube | 10fr | 120cm |
| S-NJFT-12120 | Naso Jejunal Feeding Tube | 12fr | 120cm |
| S-NJFT-14120 | Naso Jejunal Feeding Tube | 14fr | 120cm |
| S-NJFT-16120 | Naso Jejunal Feeding Tube | 16fr | 120cm |
| S-NJFT-10155 | Naso Jejunal Feeding Tube | 10fr | 155cm |
| S-NJFT-12155 | Naso Jejunal Feeding Tube | 12fr | 155cm |
| S-NJFT-14155 | Naso Jejunal Feeding Tube | 14fr | 155cm |
| S-NJFT-16155 | Naso Jejunal Feeding Tube | 16fr | 155cm |