

THE LEADER IN ENTERAL FEEDING INTRODUCES THE NEW UPGRADED MIC* PEG KIT

- New Kit & Components Streamline Procedure Placement¹
- Unique Bumper Design Helps Avoid Accidental Patient Dislodgements²



THE CHALLENGES OF DISLODGEMENT:


- According to a recent study, PEG pullouts cost more than \$20,000 in additional health care charges associated with these events⁶
- PEG pullouts often require⁶:
 - emergency department visit
 - level 3 surgical consultation
 - replacement gastrostomy tube
 - radiographic confirmation of tube positioning

THE SOLUTION:

- The MIC* PEG tube beats the competition by providing:
 - a secure fit that requires more force for removal than PEG tubes offered by Bard® and Boston Scientific®²
 - a unique tube design for simple traction removability that requires more force to remove and may result in fewer accidental dislodgments², §


RETENTION VALUES FOR PEG DEVICES

20 Fr PEG's conditioned with 1 hour 37 °C water bath.



Laboratory tests show Avanos* MIC* PEG Feeding Tubes collapse for traction removal but require more force than other tubes.²

Internal bumper features unique alternating indentations.



Bumper collapse depends on many factors. Photo is a simulation.

The bumper is designed to collapse inward upon removal.

REDUCE DISLODGEMENTS WITH OUR EXCLUSIVE TUBE DESIGN

INTUITIVE PEG KIT ORGANIZATION


- Promotes procedural efficiency, reduced procedure time & throughput¹
- Order of kit components match procedural steps of **Directions For Use**³

OPTIMIZE ESSENTIAL COMPONENTS

- Sterile Lidocaine[†] vial with rubber septum cap designed to reduce sharps risk⁴ (Safety Kit Only)
- Longer introducer needle as compared to major competitors⁵



Product Code	Description
0640-20	0640-20 PEG KIT 20F INTERNATIONALx1
0640-24	0640-24 SURG MIC ITLPEG PEG24F 1X1
0644-20	0644-20 INTERNATIONAL PULL PEG KIT 1x1
0644-24	0644-24 INTERNATIONAL PULL PEG KITx1
0645-14	0645-14 KC PEG TUBES



MIC* PEG introducer needle is 3.2 inches!
Others fall short at only 3 inches.

Comparison of lengths of PEG Safety Kit
introducer needles:

- A. BARD Ponsky[®]
- B. Boston Scientific EndoVive[®]
- C. Avanos* MIC*

Available in the U.S. only

PEG Kit



Used for percutaneous endoscopic gastrostomy placement
to provide enteral nutrition to the patient requiring nutritional support

Technical Information

Size - 20 - 22 - 24 Fr. Pull technique

Other size can be produced on request.

Availability : Sterile-for single use only.

PEG Kit CONTAINS :

1 Silicon feeding tube.

1 Loop insertion wire

1 Cable tie

1 Universal adaptor.

1 Bolus adaptor.

1 Silicon bush for locking.

Introducer needle - Big needle

11 No. scalpel

5 ml. Syringe with needle.

25 Gauge needle - Small needle.

1 Scissor.

1 Instruction manual.

Availability of snare is optional.



Ordering Information- - Mention product code x quantity required

Reference	Product Name	Size	Remarks
S-PEG-20	PEG Kit	20FR	
S-PEG-20-SN	PEG Kit	20FR	With Snare
S-PEG-22	PEG Kit	22FR	
S-PEG-22-SN	PEG Kit	22FR	With Snare
S-PEG-24	PEG Kit	24FR	
S-PEG-24-SN	PEG Kit	24FR	With Snare