

Foreign Body Holding Forcep 2 Prong, 3 Prong



Used to retrieve foreign body from G.I. tract

Technical Information :


Compatible with 2.0 & 2.8 mm. instrument channel.

Working length 170 cm. / 230 cm.

Supplied Non Sterile and It is Reusable.

Excellent and secure grip to retrieve Foreign

Body from G.I. Tract.



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Working Length	Required Channel Ø	Jaws
SFBH-223	F B Holder 2 Prong	170 cm	2.8 mm	2 Prong
SFBH-218	F B Holder 2 Prong	170 cm	2.0 mm	2 Prong
SFBH-323	F B Holder 3 Prong	170 cm	2.8 mm	3 Prong
SFBH-318	F B Holder 3 Prong	170 cm	2.0 mm	3 Prong
SFBH-423	F B Holder 4 Prong	170 cm	2.8 mm	4 Prong
SFBH-223-C	F B Holder 2 Prong	230 cm	2.8 mm	2 Prong
SFBH-218-C	F B Holder 2 Prong	230 cm	2.0 mm	2 Prong
SFBH-323-C	F B Holder 3 Prong	230 cm	2.8 mm	3 Prong
SFBH-423-C	F B Holder 4 Prong	230 cm	2.8 mm	4 Prong

“Shaili “ Range of Foreign Body Holder : 2 prong , 3 Prong , Retrieval Basket, Snare, Shark Tooth Forcep, Alligator Forcep, Monofilament Memory Basket, Loop Basket.