

Loop Basket



Used for endoscopic retrieval of foreign bodies from the G.I. tract
and to retrieve polyp through colono scope

Technical Information

Shape : Oval

With flush port.


Availability : Non sterile and for single use only.

Loop Diameter $\geq$ 61 mm could be with 2 hand operation handle

One Hand Operation



Two Hand Operation



Ordering Information- - Mention product code x quantity required

Reference	Product Name	Required Channel Ø	Loop Size	Length
SFBH-L-30	Loop Basket	2.8mm	30mm	170cm
SFBH-L-40	Loop Basket	2.8mm	40mm	170cm
SFBH-L-60	Loop Basket	2.8mm	60mm	170cm
SFBH-L-90	Loop Basket	2.8mm	90mm	170cm
SFBH-L-C-30	Loop Basket	2.8mm	30mm	230cm
SFBH-L-C-40	Loop Basket	2.8mm	40mm	230cm
SFBH-L-C-60	Loop Basket	2.8mm	60mm	230cm
SFBH-L-C-90	Loop Basket	2.8mm	90mm	230cm
SFBH-LB-DBE-60	Loop Basket	2.8mm	60mm	280cm