

Polypectomy Snare



Used for endoscopic excision of polyps using monopolar electrosurgical unit and retrieval of foreign bodies from G.I. tract

Technical Information

Availability - Sterile for Single use only.

With flush port to clean wire during procedure.

Outer sheath - PTFE.

Loop Diameter - 6 - 10 - 13 - 15 - 20 - 25 - 30 - 35 - 40 - 45 - 50 - 55 mm

Loop Diameter for Duckbill Snare is 25 mm only

Shape - Hexagonal / Oval / Crescent / Duckbill available in Multifilament

Monofilament available in only Oval shape

Individually packed .

Rotatable Snare is available



Ordering Information- Mention the product code x Loop Diameter x quantity required.

Example : Polypectomy Snare Oval x 20mm loop x Rotatable =SPSO-20-R

Reference	Product Name	Length	Shape	Required Channel Ø
SPSD - 25	Polypectomy Snare	220 cm	Duckbill	2.8 mm
SPSO - Loop Ø-R	Polypectomy Snare Rotatable	220 cm	Oval	2.8 mm
SPSH - Loop Ø-R	Polypectomy Snare Rotatable	220 cm	Hexagonal	2.8 mm
SPSC - Loop Ø-R	Polypectomy Snare Rotatable	220 cm	Crescent	2.8 mm
SPSM - Loop Ø-R	Polypectomy Snare Rotatable	220 cm	Monofilament	2.8 mm
SPSO - Loop Ø	Polypectomy Snare	220 cm	Oval	2.8 mm
SPSH - Loop Ø	Polypectomy Snare	220 cm	Hexagonal	2.8 mm
SPSC - Loop Ø	Polypectomy Snare	220 cm	Crescent	2.8 mm
SPSM - Loop Ø	Polypectomy Snare	220 cm	Monofilament	2.8 mm
SPSO -30-DBE	Polypectomy Snare	280 cm	Oval	2.8 mm
SPSO -15-SPY	Polypectomy Snare for Spyglass	280 cm	Oval	1.2 mm

Tailormade device can be produced provided by pre approved drawing/specification