

Cleaning Brush

Used to clean endoscope channel



Technical Information

Type - Non sterile single use / Reusable.

Shaft - Plastic coated wire/ Metal coil.

Semi - Soft and channel friendly bristles for Without invasive cleaning of endoscope channel.

Tip - Atraumatic.

Length - 120 cm / 180 cm /230 cm.



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Required Channel Ø	Working Length
S - CB - DE	Cleaning Brush Double Ended	2.8 mm.	230 cm
S - CB - DE - BS	Cleaning Brush Double Ended Big-Small	2.8 mm.	230 cm
S - CB - P	Port Cleaning Brush	-	-
S - CB - B	Body Cleaning Brush	-	-
S - CB - RG	Cleaning Brush Reusable Gastro	2.8 mm.	180 cm
S - CB - RC	Cleaning Brush Reusable Colono	2.8 mm.	230 cm
S - CB - D	Cleaning Brush for Dilator	1.2 mm.	120 cm