

Biopsy Forcep



Used endoscopically to obtain mucosal tissue samples from gastrointestinal tract for microscopic - histopathological examination

Technical Information


Type - With Spike / Without Spike.

Availability - Sterile, For Single use only.

Compatible with 2.8 mm instrument channel.

Length - 170cm / 210cm & 270cm

Individually packed.



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Type	Required Channel Ø	Working Length
S- BF- WS-17	Biopsy Forcep	With Spike	2.8 mm	170 cm
S- BF- WS-21	Biopsy Forcep	With Spike	2.8 mm	210 cm
S- BF- 17	Biopsy Forcep	Without Spike	2.8 mm	170 cm
S- BF- 21	Biopsy Forcep	Without Spike	2.8 mm	210 cm
S- BF- 27	Biopsy Forcep For Spy	Without Spike	1.2 mm	270 cm
S- BF- 18	Biopsy Forcep	Without Spike	2.0 mm	170 cm