

Hot Biopsy Forcep



Used endoscopically to obtain mucosal tissue samples from gastrointestinal tract for microscopic - histopathological examination

Monopolar electro surgical unit can be used to cauterize the tissue when bleeding occurs during tissue sampling

Technical Information

Without Spike

Availability- Sterile, for single use

Compatible with 2.8 instrument channel

Working Length- 170, 210cm

Active cord available separately.

Individually packed



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Type	Required Channel Ø	Working Length
S- BF-H-17	Hot Biopsy Forcep	Without Spike	2.8 mm	170 cm
S- BF-H-21	Hot Biopsy Forcep	Without Spike	2.8 mm	210 cm
S- BF-H-WS-17	Hot Biopsy Forcep	With Spike	2.8 mm	170 cm
S- BF-H-WS-21	Hot Biopsy Forcep	With Spike	2.8 mm	210 cm