

Sphincterotome

Used for endoscopic cannulation of biliary ductal system and sphincterotomy



Technical Information

Triple Lumen with Precurved tapered tip.

Tip diameter 5.5fr.

Working length 200 cm.

Wire Guide Compatibility 0.035"


Cutting wire with or without insulation as required

Protective insulation sheath, unique protective feature for endoscope,

Insulated wire length - for 20mm Loop - 4 mm, 25mm Loop - 8mm, 30 mm Loop - 10mm

Supplied Sterile for single use.

Active Cord available separately.



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Tip length	Cutting Loop	Remarks
SSPH-SN-20	Sphincterotome Short Nose	5.0mm	20mm	With Insulation
SSPH-SN-25	Sphincterotome Short Nose	5.0mm	25mm	With Insulation
SSPH-SN-30	Sphincterotome Short Nose	5.0mm	30mm	With Insulation
SSPH-SN-NI-20	Sphincterotome Short Nose	5.0mm	20mm	Without Insulation
SSPH-SN-NI-25	Sphincterotome Short Nose	5.0mm	25mm	Without Insulation
SSPH-SN-NI-30	Sphincterotome Short Nose	5.0mm	30mm	Without Insulation
SSPH-LN-20	Sphincterotome Long Nose	8.0mm	20mm	With Insulation
SSPH-LN-25	Sphincterotome Long Nose	8.0mm	25mm	With Insulation
SSPH-LN-30	Sphincterotome Long Nose	8.0mm	30mm	With Insulation
SSPH-LN-NI-20	Sphincterotome Long Nose	8.0mm	20mm	Without Insulation
SSPH-LN-NI-25	Sphincterotome Long Nose	8.0mm	25mm	Without Insulation
SSPH-LN-NI-30	Sphincterotome Long Nose	8.0mm	30mm	Without Insulation
SSPH-NN-20	Sphincterotome No Nose	0.0mm	20mm	With Insulation
SSPH-NN-25	Sphincterotome No Nose	0.0mm	25mm	With Insulation
SSPH-NN-30	Sphincterotome No Nose	0.0mm	30mm	With Insulation
SSPH-NN-NI-20	Sphincterotome No Nose	0.0mm	20mm	Without Insulation
SSPH-NN-NI-25	Sphincterotome No Nose	0.0mm	25mm	Without Insulation
SSPH-NN-NI-30	Sphincterotome No Nose	0.0mm	30mm	Without Insulation
SSPH-RC-NI-20	Sphincterotome Reverse Cut	2.0mm	20mm	Without Insulation
SSPH-RC-NI-25	Sphincterotome Reverse Cut	2.0mm	25mm	Without Insulation
SSPH-RC-NI-30	Sphincterotome Reverse Cut	2.0mm	30mm	Without Insulation
SSPH-SN-20-DBE	Sphincterotome Short Nose 280 cm.long	5.0mm	20mm	With Insulation
SSPH-SN-25-DBE	Sphincterotome Short Nose 280 cm.long	5.0mm	25mm	With Insulation
SSPH-SN-30-DBE	Sphincterotome Short Nose 280 cm long	5.0mm	30mm	With Insulation
SSPH-LN-20-DBE	Sphincterotome Long Nose 280 cm.long	8.0mm	20mm	With Insulation
SSPH-LN-25-DBE	Sphincterotome Long Nose 280 cm.long	8.0mm	25mm	With Insulation
SSPH-LN-30-DBE	Sphincterotome Long Nose 280 cm.long	8.0mm	30mm	With Insulation

Tailormade device can be produced provided by pre approved drawing/specification