

Needle Knife

Used for precut papillotomy



Technical Information

Tapered tip.

Individually packed.

Supplied Sterile and it is for single use.

Active Cord available separately.

Wire Guided Needle Knife accept 0.035" Guide Wire

Needle Knife



Wire Guided Needle Knife



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Length	Needle Length	Outer Ø
SNK-100	Needle Knife	200 cm	5 mm	1.7 mm
SNK-100-WG	Wire Guided Needle Knife	200 cm	5 mm	2.45 mm
SNK-100-DBE	Needle Knife Double Balloon Enteroscope	280 cm	5 mm	1.7 mm
SNK-100-WG-DBE	Needle Knife Double Balloon Enteroscope	280cm	5 mm	2.45 mm