


GASTRIC FEEDING TUBES

ADDITIONAL FEATURES OF THE AVANOS* G FEEDING TUBE INCLUDES:

- High clarity medical grade silicone construction
- Clearly marked feeding, medication and balloon ports
- Ventilated external retention ring
- Radiopaque stripe to aid in catheter visualization
- Graduated centimetre markings
- Color-coded, Luer-Lok balloon port



Product Code	Description
0100-12LV	0100-12LV GASTROSTOMY TUBE 12F 1x1
0100-14LV	0100-14LV GASTROSTOMY TUBE 14F 1x1
0100-16LV	0100-16LV GASTROSTOMY TUBE 16F 1x1
0100-18	0100-18 GASTROSTOMY TUBE 18F 1x1
0100-20	0100-20 GASTROSTOMY TUBE 20F 1x1
0100-22	0100-22 GASTROSTOMY TUBE 22F 1x1
0112-14 LV	0112-14LV GASTROSTOMY FEED TUBE 14F 1X1
0112-18	0112-18 GASTROSTOMY FEED TUBE 18FR 1x1
0112-20	0112-20 GASTROSTOMY FEED TUBE 20FR 1x1
0112-22	0112-22 GASTROSTOMY FEED TUBE 22FR 1x1
0112-24	0112-24 GASTROSTOMY FEED TUBE 24FR 1x1
0112-16 LV	0112-16LV GASTROSTOMY FEED TUBE 16F 1X1

1. Data on file. DVC-2013/-35 - Design Inputs.
2. Data on file. GL-DSR00143 / 1 - Retention Values for Competitor PEG Tubes Report.
3. Data on file. MIC* Pull PEG Kit Directions for Use, MIC* Push PEG Kit Directions for Use, MIC* Safety Pull PEG Kit Directions for Use & MIC* Safety Push PEG Kit Directions for Use.
4. Carraretto AR, et al. Glass Ampoules: Risks and Benefits. Rev Bras Anestesiol, 2011; 61: 4: 513-521 & Data on file. Avanos Medical.
5. Data on file. Table of Competitive Needle Lengths from PEG Placement Kit - Highlights and Market Comparison.
6. Rosenberger LH, et al. Late accidental dislodgement of a percutaneous endoscopic gastrostomy tube: an underestimated burden on patients and the health care system. Surg Endosc. 2011 October; 25(10): 3307-3311

For more information please visit: avanosmedicaldevices.com
Call 1-844-4AVANOS (1-844-428-2667) in the United States and Canada.

Replacement PEG Tube



Used as a replacement gastrostomy tube through an established gastrostomy treatment

Technical Information

Size - 20 - 22 - 24 Fr.

Availability : Sterile-for single use only.

- 1 pc Gastrostomy Replacement Tube
- 1 pc 12 ml syringe
- 2 pc 10 ml Normal saline



Ordering Information- - Mention product code x quantity required

Reference	Product Name	Size
S - GRT - 20	Gastrostomy Replacement Tube	20 Fr.
S - GRT - 22	Gastrostomy Replacement Tube	22 Fr.
S - GRT - 24	Gastrostomy Replacement Tube	24 Fr.