

Endo Scissor

Used to cut and coagulate during
endoscopic mucosal resection



Technical Information

Working length 160/230cm

Sheath diameter 2.5 mm

Availability - Sterile for single use only.

Rotatable scissor is available with diathermy connection

Active Cord available separately



Ordering Information : Mention the product code x quantity required.

Reference	Product Name	Length	Sheath Diameter
S-ES-160-H	Endo Scissor Hot	160cm	2.5mm
S-ES-160-R-H	Endo Scissor Hot Rotatable	160cm	2.5mm
S-ES-230-H	Endo Scissor Hot	230cm	2.5 mm
S-ES-230-R-H	Endo Scissor Hot Rotatable	230cm	2.5 mm