

Retrieval Basket



Used to retrieve foreign bodies and C.B.D. stones endoscopically

Technical Information :

Compatible with instrument channel Ø 2.8 mm.

Basket diameter - 15 mm. / 20 mm. / 25 mm. / 30 mm. / 35 mm. / 40 mm. / 45 mm. / 50 mm.


Multifilament, Hexagonal shaped.

Flush port for contrast.

Atraumatic basket tip for ease of insertion and procedure.

Compatible with Lithotripsy.

Sterile for single use only.



Ordering Information- Strand diameter is 0.41mm , Working length 230 cm.

Reference	Product Name	Basket Ø	Required Channel Ø	Working Length
SRB-4152	Retrieval Basket	15 mm	2.8 mm	230 cm
SRB-4202	Retrieval Basket	20 mm	2.8 mm	230 cm
SRB-4252	Retrieval Basket	25 mm	2.8 mm	230 cm
SRB-4302	Retrieval Basket	30 mm	2.8 mm	230 cm
SRB-4352	Retrieval Basket	35 mm	2.8 mm	230 cm
SRB-4402	Retrieval Basket	40 mm	2.8 mm	230 cm
SRB-4452	Retrieval Basket	45 mm	2.8 mm	230 cm
SRB-4502	Retrieval Basket	50 mm	2.8 mm	230 cm
SRB-4302-DBE	Retrieval Basket	30 mm	2.8 mm	280 cm

Tailormade device can be produced provided by pre approved drawing/specification