

# Wire Guided Retrieval Basket



Used to retrieve foreign bodies and C.B.D. stones endoscopically

## Technical Information :

Working length 230 cm.

Separate lumen for wire.

Compatible with instrument channel 3.2 mm.

Basket diameter -15mm. / 20 mm. / 25mm. / 30 mm./35 mm.

40 mm./ 45 mm./ 50 mm.

Multifilament, Hexagonal shaped.

Flush port for contrast.

Atraumatic basket tip for ease of insertion and procedure.

**Compatible with Lithotripsy.**

Sterile for single use only.



**Ordering Information- Strand diameter is 0.41mm, Working length 230 cm.**

| Reference   | Product Name                 | Basket Ø | Sheath Ø | Wire Guide Compatibility |
|-------------|------------------------------|----------|----------|--------------------------|
| SRB-WG-4153 | Wire Guided Retrieval Basket | 15 mm    | 3.0 mm   | 0.035"                   |
| SRB-WG-4203 | Wire Guided Retrieval Basket | 20 mm    | 3.0 mm   | 0.035"                   |
| SRB-WG-4253 | Wire Guided Retrieval Basket | 25 mm    | 3.0 mm   | 0.035"                   |
| SRB-WG-4303 | Wire Guided Retrieval Basket | 30 mm    | 3.0 mm   | 0.035"                   |
| SRB-WG-4353 | Wire Guided Retrieval Basket | 35 mm    | 3.0 mm   | 0.035"                   |
| SRB-WG-4403 | Wire Guided Retrieval Basket | 40 mm    | 3.0 mm   | 0.035"                   |
| SRB-WG-4453 | Wire Guided Retrieval Basket | 45 mm    | 3.0 mm   | 0.035"                   |
| SRB-WG-4503 | Wire Guided Retrieval Basket | 50 mm    | 3.0 mm   | 0.035"                   |

**Tailormade device can be produced provided by pre approved drawing/specification**