

Tip Guided Retrieval Basket



Used to retrieve foreign bodies and C.B.D. stones

Technical Information :

Working length 230 cm.

Wire access point on tip of basket.

Compatible with instrument channel 2.8 mm or more.


Basket diameter -15mm. / 20 mm. / 25mm. / 30 mm. / 35 mm. / 40 mm. / 45 mm. / 50 mm.

Multifilament, Hexagonal shaped.

Flush port for contrast.

Compatible with Lithotripsy.

Sterile for single use only.



Ordering Information- Strand diameter is 0.41mm, Working length 230 cm.

Reference	Product Name	Basket Ø	Sheath Ø	Wire Guide Compatibility
SRB-WG-4152-T	Tip Guided Retrieval Basket	15 mm	2.5 mm	0.035"
SRB-WG-4202-T	Tip Guided Retrieval Basket	20 mm	2.5 mm	0.035"
SRB-WG-4252-T	Tip Guided Retrieval Basket	25 mm	2.5 mm	0.035"
SRB-WG-4302-T	Tip Guided Retrieval Basket	30 mm	2.5 mm	0.035"
SRB-WG-4352-T	Tip Guided Retrieval Basket	35 mm	2.5 mm	0.035"
SRB-WG-4402-T	Tip Guided Retrieval Basket	40 mm	2.5 mm	0.035"
SRB-WG-4452-T	Tip Guided Retrieval Basket	45 mm	2.5 mm	0.035"
SRB-WG-4502-T	Tip Guided Retrieval Basket	50 mm	2.5 mm	0.035"

Tailormade device can be produced provided by pre approved drawing/specification