

Memory Basket

Used to retrieve foreign body / stone from G.I. Tract - through the scope

Technical Information

Seven basket sizes to facilitate foreign body retrieval from G.I. Tract / stone retrieval from bileduct.


Basket Type - Oval, Hexagonal, Tipless & Helical

Retains shape for longer time.

Handle - Light weight with flush port for contrast-drug-flush.

Not recommended for Lithotripsy.

Availability - Sterile - for single use only



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Required Channel Ø	Basket Ø	Working Length
SMMB-15	Memory Basket	2.8mm	15mm	230cm
SMMB-20	Memory Basket	2.8mm	20mm	230cm
SMMB-25	Memory Basket	2.8mm	25mm	230cm
SMMB-30	Memory Basket	2.8mm	30mm	230cm
SMMB-35	Memory Basket	2.8mm	35mm	230cm
SMMB-40	Memory Basket	2.8mm	40mm	230cm
SMMB-45	Memory Basket	2.8mm	45mm	230cm
SMMB-15-T	Memory Basket Tipless	2.8mm	15mm	230cm
SMMB-20-T	Memory Basket Tipless	2.8mm	20mm	230cm
SMMB-25-T	Memory Basket Tipless	2.8mm	25mm	230cm
SMMB-30-T	Memory Basket Tipless	2.8mm	30mm	230cm
SMMB-35-T	Memory Basket Tipless	2.8mm	35mm	230cm
SMMB-40-T	Memory Basket Tipless	2.8mm	40mm	230cm
SMMB-45-T	Memory Basket Tipless	2.8mm	45mm	230cm

Tailormade device can be produced provided by pre approved drawing/specification