

Wire Guided Memory Basket

Used to retrieve foreign body / stone from G.I. Tract - through the scope

Technical Information

Three basket sizes to facilitate stone retrieval from biliary tract.

Basket shapes - Oval, Hexagonal.

Retains shape for longer time.

Handle - Light weight with flush port for contrast.

Centre port for guide wire

Wire guide compatibility 0.035"

Not recommended for Lithotripsy.

Availability - Sterile - for single use only



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Required Channel Ø	Basket Ø	Working Length
SMMB-WG-30	Wire Guided Memory Basket	2.8mm	30mm	230cm
SMMB-WG-35	Wire Guided Memory Basket	2.8mm	35mm	230cm
SMMB-WG-40	Wire Guided Memory Basket	2.8mm	40mm	230cm

Tailormade device can be produced provided by pre approved drawing/specification