

Litho Basket-TTS



Used for purpose of C.B.D. stone retrieval/ lithotripsy of stone in bile duct through the scope

Technical Specification :

Working length 230cm

Extension wire length 205cm


Loop Diameter 25/30/35/40/45/50/55mm

Sheath Diameter 2.5mm

Recommended for manual lithotripsy.

Supplied sterile and it is for single use

Detachable Handle.



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Loop Ø	Sheath Ø	Working Length
SLTB-25	Lithotripsy Basket TTS	25mm	2.5mm	230cm
SLTB-30	Lithotripsy Basket TTS	30mm	2.5mm	230cm
SLTB-35	Lithotripsy Basket TTS	35mm	2.5mm	230cm
SLTB-40	Lithotripsy Basket TTS	40mm	2.5mm	230cm
SLTB-45	Lithotripsy Basket TTS	45mm	2.5mm	230cm
SLTB-50	Lithotripsy Basket TTS	50mm	2.5mm	230cm
SLTB-55	Lithotripsy Basket TTS	55mm	2.5mm	230cm