

Manual Lithotripter

Used for C.B.D. stone lithotripsy




Technical Information

Handle - Light weight with unidirectional movement.

Sheath tip - Inner cutting edges to perform lithotripsy. Outer tip rounded atraumatic.

Secure and firm fitting of S.S. Sheath with handle to facilitate retrieval - lithotripsy procedure.

Availability - Non Sterile Reusable.



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Spring Sheath Ø	Spring Sheath Length	Required Channel Ø
SML-SHEATH-16590	Manual Lithotripter Spring Sheath	16.5FR	90cm	NA
SML-SHEATH-1190	Manual Lithotripter Spring Sheath	11 FR	90cm	NA
SML-SHEATH-11180-TTS	Manual Lithotripter TTS Spring Sheath	11 FR	180cm	4.2mm
SML-SHEATH-75180-TTS	Manual Lithotripter TTS Spring Sheath	7.5 FR	180cm	2.8mm
SML-SHEATH-9180-TTS	Manual Lithotripter TTS Spring Sheath	9 FR	180cm	3.2mm
SML-H	Manual Lithotripter Handle only			