

Nitinol Black Wire



Used to gain access through difficult biliary and pancreatic strictures.
To facilitate exchange of compatible wire guided accessories to the
desired biliary or pancreatic ducts

Technical Information

Size Diameter - 0.025 inch / 0.032 inch / 0.035 inch. / 0.038 inch.


Availability- Individually packed. Sterile For single use only.

Length - 260 cm / 450 cm

Tip - Straight / Angle / J

Hydrophilic Coated

Straight Angle J



Ordering Information : Mention the product code x quantity required.

Reference	Product Name	Reference	Product Name
S-NBW-025-S X LENGTH	Nitinol Black Wire 0.025 Straight	S-NBW-035-S X LENGTH	Nitinol Black Wire 0.035 Straight
S-NBW-025-A X LENGTH	Nitinol Black Wire 0.025 Angulated	S-NBW-035-A X LENGTH	Nitinol Black Wire 0.035 Angulated
S-NBW-025-J X LENGTH	Nitinol Black Wire 0.025 J Tip	S-NBW-035-J X LENGTH	Nitinol Black Wire 0.035 J Tip
S-NBW-032-S X LENGTH	Nitinol Black Wire 0.032 Straight	S-NBW-038-S X LENGTH	Nitinol Black Wire 0.038 Straight
S-NBW-032-A X LENGTH	Nitinol Black Wire 0.032 Angulated	S-NBW-038-A X LENGTH	Nitinol Black Wire 0.038 Angulated
S-NBW-032-J X LENGTH	Nitinol Black Wire 0.032 J Tip	S-NBW-038-J X LENGTH	Nitinol Black Wire 0.038 J Tip