

Wire Guided Stone Extraction Balloon



Used for endoscopic extraction of biliary stones

Technical Information

Technical Information-

Triple lumen feature provides wire guide access to facilitate stone extraction, from desired location during procedure.

Catheter Diameter - 7 Fr .


Compatible with 0.035 inch wire guide.

Availability - Individually packed. Sterile for single use only.

Working length - 200 cm.

Balloon diameter - 16 mm.

Balloon length - 16 mm.



Reference	Catheter Diameter FR	Functional Catheter Length	Diameter of Inflated Balloon	Balloon Max Air Inflation	Balloon Max Fluid Injection	Recommended Guiding Wire (Inch)	Syringe CC	Tip Diameter FR
SSEC - 7	7	200 cm	16 mm	3.5 cc	3.0 ml	0.035"	5	7