

Achalasia Balloon Dilator



Intended for dilation of the lower esophageal sphincter
in patient with achalasia

Technical Information :

Inflated Balloon Diameter - 25mm, 30mm, 35 mm

Balloon length : 80 mm.

Shaft Diameter : 5mm.

Working Length : 90 cm.

Wire Guide Compatibility : 0.038"


Burst pressure : 3 ATM

High Radial Vector Force.

Atraumatic Tip.

Radioopaque marker inside the balloon.

Sterile for Single use only.



Ordering Information - Mention product code x quantity required

Reference	Product Name	Balloon Ø
SABD - 25	Achalasia Balloon Dilator	25 mm
SABD - 30	Achalasia Balloon Dilator	30 mm
SABD - 35	Achalasia Balloon Dilator	35 mm