

OmniView MULTIPLE BAND LIGATOR KIT



A variceal ligation system used for endoscopic ligation of esophageal varices, at or above the gastroesophageal junction.

PACKAGING: 01 piece per box, consisting of:

- Ligator cartridge with ultra-clear, wide view polycarbonate barrel preloaded with high quality elastomeric bands. **CAUTION:** The elastomeric bands are made of natural rubber latex.
- Deployment handle, with an ultra smooth 'one-way' movement, for reliable band deployment.
- Loading catheter
- Irrigation needle

The product is available in various sizes to fit endoscopes ranging from 8.6 mm to 14 mm diameter.

MBLS-4	Set consisting of 1 cartridge with 4 bands for endoscope diameter 9.5 ~ 13.0 mm
MBLS-6	Set consisting of 1 cartridge with 6 bands for endoscope diameter 9.5 ~ 13.0 mm
MBLS-4-XL	Set consisting of 1 cartridge with 4 bands for endoscope diameter 11.0 ~ 14.0 mm
MBLS-6-XL	Set consisting of 1 cartridge with 6 bands for endoscope diameter 11.0 ~ 14.0 mm
MBLS-4-XS	Set consisting of 1 cartridge with 4 bands for endoscope diameter 8.6 ~ 9.2 mm
MBLS-6-XS	Set consisting of 1 cartridge with 6 bands for endoscope diameter 8.6 ~ 9.2 mm
MBLS-6-OV	OptiVue Set consisting of 1 cartridge with 6 bands for endoscope diameter 9.5 ~ 11.0 mm
MBLS-6-XS-OV	OptiVue Set consisting of 1 cartridge with 6 bands for endoscope diameter 8.6 ~ 9.2 mm

This product is supplied non-sterile and intended for single use, to be used by or on order of a physician only.

Marketed by

Globax Marketing & Consultants

303, Third Floor, Somdutt Chamber - II, 9, Bhikaji Cama Place, New Delhi - 110066

T: 011-49051596/26176463, Website: www.globax.in, Email: info@globax.in

Medelec
SYSTEMS PVT LTD

Multi Band Ligator Barrel with 6 loaded Bands.



Used for Endoscopic Ligation of Esophageal Varices at or above Gastro Esophageal Junction.

This device is supplied non sterile, individually packed - intended for single use only.

Technical Information :

Shaili Band Ligator requires 9.5 to 11.2 .mm endoscope diameter.


Shaili Paediatric Band Ligator requires 8.0 to 9.5 mm endoscope diameter.

Band Ligator Set Contains

Applicator x 1 pc., Barrel loaded with 6 x Bands,
Barrel loaded with 7 x Bands,
Atraumatic plastic hook with Tube X1pc.,
Flushing Needle blunt tip x1 pc.
Second last band with different colour.

Barrel Contains

Barrel loaded with 6 x Bands,
Barrel loaded with 7 x Bands,
Atraumatic plastic hook with Tube X1pc.



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Reference	Product Name
SBLS-6-F*	Multi Band Ligator Set	SBLS-6-O*	Multi Band Ligator Set
SBLS-7-F*	Multi Band Ligator Set	SBLS-7-O*	Multi Band Ligator Set
SBL-6-F*	Multi Band Ligator Barrel	SBL-6-O*	Multi Band Ligator Barrel
SBL-7-F*	Multi Band Ligator Barrel	SBL-7-O*	Multi Band Ligator Barrel
SBLS-6-P-F*	Multi Band Ligator Set Paediatric	SBLS-6-P-O*	Multi Band Ligator Set Paediatric
SBL-6-P-F*	Multi Band Ligator Barrel Paediatric	SBL-6-P-O*	Multi Band Ligator Barrel Paediatric
SBLS-7-P-F*	Multi Band Ligator Set Paediatric	SBLS-7-P-O*	Multi Band Ligator Set Paediatric
SBL-7-P-F*	Multi Band Ligator Barrel Paediatric	SBL-7-P-O*	Multi Band Ligator Barrel Paediatric

* To Indicate the Length of Thread., F = 175 cm. for Fujinon Endoscope, O = 155 cm. For Olympus / Pentax Endoscope.

Multi Band Ligator Colonic Barrel with 6 loaded Bands.



Used for endoscopic ligation of colonic varices

This device is supplied non sterile, individually packed - intended for single use only.

Technical Information :

Shaili Band Ligator Colonic requires 12.8 mm to 13.5 mm endoscope diameter.

Band Ligator Colonic Set Contains

Applicator x 1 pc., Barrel loaded with 6 x Bands,
Barrel loaded with 7 x Bands,
Atraumatic plastic hook with Tube X1pc.,
Flushing Needle blunt tip x1 pc.
Second last band with different colour.

Barrel Colonic Contains

Barrel loaded with 6 x Bands,
Barrel loaded with 7 x Bands,
Atraumatic plastic hook with Tube X1pc.



Ordering Information- Mention the product code x quantity required.

Reference	Product Name
SBLS-6-C	Multi Band Ligator Colonic Set
SBL-6-C	Multi Band Ligator Colonic Barrel
SBLS-7-C	Multi Band Ligator Colonic Set
SBL-7-C	Multi Band Ligator Colonic Barrel