

# Biliary Stent

Used to drain bile from obstructed biliary duct



## Technical Information

Diameter 5 Fr. - 7Fr. - 8.5 Fr. - 10 Fr. & 11.5 Fr. Size

Length 3 cm. to 20 cm. as per customer requirement.

Material choice- Radioopaque PTFE OR PE.

Varients


- Straight - S
- Double Pigtail - DP
- Amsterdam - AMST
- Standard - STD
- S Banded - SB
- Single Pigtail - SP
- C Shape - CS
- Tree Type (Pine tree) - T

Available - Sterile, Single use.

Wireguide Compatibility - 0.035".

TEF = PTFE material

PU = PE material



**Ordering Information- Mention the product code x required Length x QTY required.**

**Example : Biliary Stent X PE X 7fr X Standard X10cm length = SBS-PU-7-STD-10**

| Specifications  | Stent 5 Fr. Reference | Biliary Stent 7 Fr. Reference | Biliary Stent 8.5 Fr. Reference | Biliary stent 10Fr. Reference | Biliary stent 11.5 Fr. Reference |
|-----------------|-----------------------|-------------------------------|---------------------------------|-------------------------------|----------------------------------|
| Material        |                       |                               |                                 |                               |                                  |
| PTFE            | SBS-TEF-5             | SBS-TEF-7                     | SBS-TEF-8.5                     | SBS-TEF-10                    | SBS-TEF-11.5                     |
| PE              | SBS-PU-5              | SBS-PU-7                      | SBS-PU-8.5                      | SBS-PU-10                     | SBS-PU-11.5                      |
| Variant :       |                       |                               |                                 |                               |                                  |
| Straight        | SBS-TEF-5-S           | SBS-TEF-7- S                  | SBS-TEF-8.5-S                   | SBS-TEF-10-S                  | SBS-TEF-11.5-S                   |
|                 | SBS-PU-5-S            | SBS-PU-7-S                    | SBS-PU-8.5-S                    | SBS-PU-10-S                   | SBS-PU-11.5-S                    |
| Double Pig Tail | SBS-PU-5-DP           | SBS-PU-7-DP                   | SBS-PU-8.5-DP                   | SBS-PU-10-DP                  | SBS-PU-11.5-DP                   |
| Standard        | SBS-TEF-5-STD         | SBS-TEF-7-STD                 | SBS-TEF-8.5-STD                 | SBS-TEF-10-STD                | SBS-TEF-11.5-STD                 |
|                 | SBS-PU-5-STD          | SBS-PU-7-STD                  | SBS-PU-8.5-STD                  | SBS-PU-10-STD                 | SBS-PU-11.5-STD                  |
| Amsterdam       | SBS-TEF-5-AMST        | SBS-TEF-7-AMST                | SBS-TEF-8.5-AMST                | SBS-TEF-10-AMST               | SBS-TEF-11.5-AMST                |
|                 | SBS-PU-5-AMST         | SBS-PU-7-AMST                 | SBS-PU-8.5-AMST                 | SBS-PU-10-AMST                | SBS-PU-11.5-AMST                 |
| S Banded        | SBS-PU-5-SB           | SBS-PU-7-SB                   | SBS-PU-8.5-SB                   | SBS-PU-10-SB                  | SBS-PU-11.5-SB                   |
| Single Pigtail  | SBS-PU-5-SP           | SBS-PU-7-SP                   | SBS-PU-8.5-SP                   | SBS-PU-10-SP                  | SBS-PU-11.5-SP                   |
| C Shape         | SBS-PU-5-CS           | SBS-PU-7-CS                   | SBS-PU-8.5-CS                   | SBS-PU-10-CS                  | SBS-PU-11.5-CS                   |
|                 | SBS-TEF-5-CS          | SBS-TEF-7-CS                  | SBS-TEF-8.5-CS                  | SBS-TEF-10-CS                 | SBS-TEF-11.5-CS                  |
| Tree Type       |                       |                               | SBS-TEF-8.5-T                   | SBS-TEF-10-T                  | SBS-TEF-11.5-T                   |