

Biliary Stent Set

Used to drain bile from obstructed biliary duct



Technical Information

Biliary Stent Set Contain - Biliary Stent & Pusher of 5 Fr. & 7Fr.

Biliary Stent Set 8.5 Fr. - 10 Fr. & 11.5 Fr. Contains - Biliary Stent, Pusher & Guiding Catheter


Length of Guiding Catheter 220 cm / 320 cm

Available - Sterile, Single use.

Wireguide Compatibility - 0.035".

TEF = PTFE material

PU = PE material



Ordering Information- Mention the product code x required Length x QTY required.

Example : Biliary Stent set X PE X 10 fr X Double Pigtail X 10cm long = SBSS-PU-10-DP-10

Specifications	Stent 5 Fr. Reference	Biliary Stent 7 Fr. Reference	Biliary Stent 8.5 Fr. Reference	Biliary stent 10Fr. Reference	Biliary stent 11.5 Fr. Reference
Material					
PTFE	SBSS-TEF-5	SBSS-TEF-7	SBSS-TEF-8.5	SBSS-TEF-10	SBSS-TEF-11.5
PE	SBSS-PU-5	SBSS-PU-7	SBSS-PU-8.5	SBSS-PU-10	SBSS-PU-11.5
Variant :					
Straight	SBSS-TEF-5-S	SBSS-TEF-7- S	SBSS-TEF-8.5-S	SBSS-TEF-10-S	SBSS-TEF-11.5-S
	SBSS-PU-5-S	SBSS-PU-7-S	SBSS-PU-8.5-S	SBSS-PU-10-S	SBSS-PU-11.5-S
Double Pig Tail	SBSS-PU-5-DP	SBSS-PU-7-DP	SBSS-PU-8.5-DP	SBSS-PU-10-DP	SBSS-PU-11.5-DP
Standard	SBSS-TEF-5-STD	SBSS-TEF-7-STD	SBSS-TEF-8.5-STD	SBSS-TEF-10-STD	SBSS-TEF-11.5-STD
	SBSS-PU-5-STD	SBSS-PU-7-STD	SBSS-PU-8.5-STD	SBSS-PU-10-STD	SBSS-PU-11.5-STD
Amsterdam	SBSS-TEF-5-AMST	SBSS-TEF-7-AMST	SBSS-TEF-8.5-AMST	SBSS-TEF-10-AMST	SBSS-TEF-11.5-AMST
	SBSS-PU-5-AMST	SBSS-PU-7-AMST	SBSS-PU-8.5-AMST	SBSS-PU-10-AMST	SBSS-PU-11.5-AMST
S Banded	SBSS-PU-5-SB	SBSS-PU-7-SB	SBSS-PU-8.5-SB	SBSS-PU-10-SB	SBSS-PU-11.5-SB
Single Pigtail	SBSS-PU-5-SP	SBSS-PU-7-SP	SBSS-PU-8.5-SP	SBSS-PU-10-SP	SBSS-PU-11.5-SP
C Shape	SBSS-PU-5-CS	SBSS-PU-7-CS	SBSS-PU-8.5-CS	SBSS-PU-10-CS	SBSS-PU-11.5-CS
	SBSS-TEF-5-CS	SBSS-TEF-7-CS	SBSS-TEF-8.5-CS	SBSS-TEF-10-CS	SBSS-TEF-11.5-CS
Tree Type			SBSS-TEF-8.5-T	SBSS-TEF-10-T	SBSS-TEF-11.5-T