

Pancreatic Stent

Used to drain the pancreatic juice from obstructed pancreatic duct.

Technical Information

Diameter 5 Fr. - 7Fr. - 8.5 Fr. - 10 Fr. & 11.5 Fr. Size

Length 3 cm. to 20 cm. as per customer requirement.

Material choice- Radioopaque PTFE OR PE.

Varients

Straight - S

Single Pigtail - SP


C Shape - CS

Available - Sterile, Single use.

Wireguide Compatibility - 0.035".

TEF = PTFE material

PU = PE material



Ordering Information- Mention the product code x required Length x QTY required.

Example : Pancreatic Stent X PE X 7fr X Single Pigtail X10cm length = SPS-PU-7-SP-10

Specifications	Pancreatic Stent 5 Fr. Reference	Pancreatic Stent 7 Fr. Reference	Pancreatic Stent 8.5 Fr. Reference	Pancreatic Stent 10Fr. Reference	Pancreatic Stent 11.5 Fr. Reference
Material					
PTFE	SPS-TEF-5	SPS-TEF-7	SPS-TEF-8.5	SPS-TEF-10	SPS-TEF-11.5
PE	SPS-PU-5	SPS-PU-7	SPS-PU-8.5	SPS-PU-10	SPS-PU-11.5
Variant :					
Straight	SPS-TEF-5-S	SPS-TEF- 7- S	SPS-TEF-8.5-S	SPS-TEF-10-S	SPS-TEF-11.5-S
Single Pigtail	SPS-PU-5-SP	SPS-PU- 7 - SP	SPS-PU-8.5-SP	SPS-PU-10 -SP	SPS-PU-11.5-SP
C Shape	SPS-PU-5-CS	SPS-PU- 7 - CS	SPS-PU-8.5-CS	SPS-PU-10 -CS	SPS-PU-11.5-CS
	SPS-TEF-5-CS	SPS -TEF- 7- CS	SPS-TEF-8.5-CS	SPS-TEF-10-CS	SPS-TEF-11.5-CS