

# Pancreatic Stent Set

Used to drain the pancreatic juice from obstructed pancreatic duct.



## Technical Information

Pancreatic Stent Set Contain - Pancreatic Stent & Pusher of 5 Fr. & 7Fr.

Pancreatic Stent Set 8.5 Fr. -10 Fr. & 11.5 Fr. Contains - Pancreatic Stent, Pusher & Guiding Catheter


Length of Guiding Catheter 220 cm / 320 cm

Available - Sterile, Single use.

**Wireguide Compatibility - 0.035".**

TEF = PTFE material

PU = PE material



**Ordering Information- Mention the product code x required Length x QTY required.**

**Example : Pancreatic Stent set X PE X 10 fr X Single Pigtail X 10cm long = SPSS-PU-10-SP-10**

| Specifications | Pancreatic Stent<br>5 Fr. Reference | Pancreatic Stent<br>7 Fr. Reference | Pancreatic Stent<br>8.5 Fr. Reference | Pancreatic Stent<br>10Fr. Reference | Pancreatic Stent<br>11.5 Fr. Reference |
|----------------|-------------------------------------|-------------------------------------|---------------------------------------|-------------------------------------|----------------------------------------|
| Material       |                                     |                                     |                                       |                                     |                                        |
| PTFE           | SPSS-TEF-5                          | SPSS-TEF-7                          | SPSS-TEF-8.5                          | SPSS-TEF-10                         | SPSS-TEF-11.5                          |
| PE             | SPSS-PU-5                           | SPSS-PU-7                           | SPSS-PU-8.5                           | SPSS-PU-10                          | SPSS-PU-11.5                           |
| Variant :      |                                     |                                     |                                       |                                     |                                        |
| Straight       | SPSS-TEF-5-S                        | SPSS-TEF- 7- S                      | SPSS-TEF-8.5-S                        | SPSS-TEF-10-S                       | SPSS-TEF-11.5-S                        |
| Single Pigtail | SPSS-PU-5-SP                        | SPSS-PU- 7 - SP                     | SPSS-PU-8.5-SP                        | SPSS-PU-10 -SP                      | SPSS-PU-11.5-SP                        |
| C Shape        | SPSS-PU-5-CS                        | SPSS-PU- 7 - CS                     | SPSS-PU-8.5-CS                        | SPSS-PU-10 -CS                      | SPSS-PU-11.5-CS                        |
|                | SPSS-TEF-5-CS                       | SPSS -TEF- 7- CS                    | SPSS-TEF-8.5-CS                       | SPSS-TEF-10-CS                      | SPSS-TEF-11.5-CS                       |