

Single Action Stent Introducer System



Used for endoscopic biliary stent placement
to reduce time. Easier stent placement

Technical Information

Size - 8.5fr., 10fr. & 11.5fr. Single action stent deployment.

Radioopaque Marker Over Guiding Catheter.

Low Friction Movement.


Compatible With 0.035 Inch Ø Wire Guide.

Supply Without Guide Wire and Stent .

Availability - Sterile, For Single Use Only.

Pusher length - 180cm

Guiding catheter length - 220cm



Ordering Information- - Mention product code x quantity required

Reference	Product Name	Size	Remark
SSIS - 8.5	Single Action Stent Introducer System	8.5 Fr.	Without Guide Wire and Stent.
SSIS - 10	Single Action Stent Introducer System	10 Fr.	Without Guide Wire and Stent.
SSIS - 11.5	Single Action Stent Introducer System	11.5 Fr.	Without Guide Wire and Stent.