

Single Action Stent Introducer System Set



Used for endoscopic biliary stent placement
to reduce time. Easier stent placement

Technical Information

Size - 8.5fr., 10fr. & 11.5fr. Single action stent deployment.

Radioopaque Marker Over Guiding Catheter.

Low Friction Movement.

Compatible With 0.035 Inch Ø Wire Guide.

Pusher length - 180cm

Guiding catheter length - 220cm

Supply with Biliary Stent .


Supply without Guide Wire.

Availability - Sterile, For Single Use Only.

Biliary Stent Length - 3 cm to 20 cm

TEF = PTFE material

PU = PE material



Ordering Information- Mention the product code x required Length x QTY required.

Example : single action stent introducer system set X 10 fr X Teflon X 10 fr x Straight X 10cm long = SSISS-10-TEF-10-S-10

Specifications	Single Action Stent Introducer Set 8.5 Fr.	Single Action Stent Introducer Set 10 Fr.	Single Action Stent Introducer Set 11.5 Fr.
Material			
PTFE	SSISS-8.5-TEF-8.5	SSISS-10-TEF-10	SSISS-11.5-TEF-11.5
PE	SSISS-8.5-PU-8.5	SSISS-10-PU-10	SSISS-11.5-PU-11.5
Variant :			
Straight	SSISS-8.5-TEF-8.5-S	SSISS-10-TEF-10-S	SSISS-11.5-TEF-11.5-S
	SSISS-8.5-PU-8.5-S	SSISS-10-PU-10-S	SSISS-11.5-PU-11.5-S
Double Pig Tail	SSISS-8.5-PU-8.5-DP	SSISS-10-PU-10 -DP	SSISS-11.5-PU-11.5 - DP
Standard	SSISS-8.5-TEF-8.5-STD	SSISS-10-TEF-10 -STD	SSISS-11.5-TEF-11.5-STD
	SSISS-8.5-PU-8.5-STD	SSISS-10-PU-10 -STD	SSISS-11.5-PU-11.5-STD
Amsterdam	SSISS-8.5-TEF-8.5-AMST	SSISS-10-TEF-10 -AMST	SSISS-11.5-TEF-11.5-AMST
	SSISS-8.5-PU-8.5-AMST	SSISS-10-PU-10 -AMST	SSISS-11.5-PU-11.5-AMST
Pancreatic	SSISS-8.5-TEF-8.5-PAN	SSISS-10-TEF-10-PAN	SSISS-11.5-TEF-11.5-PAN
	SSISS-8.5-PU-8.5-PAN	SSISS-10-PU-10-PAN	SSISS-11.5-PU-11.5-PAN
S Banded	SSISS-8.5-PU-8.5-SB	SSISS-10-PU-10 -SB	SSISS-11.5-PU-11.5 - SB
Single Pigtail	SSISS-8.5-PU-8.5-SP	SSISS-10-PU-10 -SP	SSISS-11.5-PU-11.5 - SP
C Shape	SSISS-8.5-PU-8.5-CS	SSISS-10-PU-10 -CS	SSISS-11.5-PU-11.5 - CS
	SSISS-8.5-TEF-8.5-CS	SSISS-10-TEF-10-CS	SSISS-11.5-TEF-11.5-CS