

Guiding Catheter

Used to place biliary stent endoscopically




Technical Information

Size - 5 fr. - 6 fr. - 7 fr.

Working Length - 45 cm, 220 cm & 320 cm

Availability - Sterile - Single use only.

Individually packed.



Ordering Information - Mention the product code x quantity required.

Reference	Product Name	Size	Working Length	Compatible Pusher
SGC-545	Guiding Catheter	5 Fr	45 cm	
SGC-645	Guiding Catheter	6 Fr	45 cm	
SGC-745	Guiding Catheter	7 Fr	45 cm	
SGC-5320	Guiding Catheter	5 Fr	320 cm	8.5 Fr
SGC-5220	Guiding Catheter	5 Fr	220 cm	8.5 Fr
SGC-6320	Guiding Catheter	6 Fr	320 cm	10 Fr
SGC-6220	Guiding Catheter	6 Fr	220 cm	10 Fr
SGC-7320	Guiding Catheter	7 Fr	320 cm	11.5 Fr
SGC-7220	Guiding Catheter	7 Fr	220 cm	11.5 Fr