

Stent Pusher

Used Endoscopically for biliary stent placement




Technical Information

Size - 5 fr. - 7 fr. - 8.5 fr. - 10fr. - 11.5fr.

Availability - Sterile - Single use only.

Individually packed.

Length : 160 cm / 180 cm / 200 cm / 230 cm / 250 cm



Ordering Information - Mention the product code x quantity required.

Reference	Product Name	Size	Required Channel Ø
SSP-5 -X Length	Stent Pusher	5 Fr	≥3.2 mm
SSP-7-X Length	Stent Pusher	7 Fr	≥3.2 mm
SSP-8.5-X Length	Stent Pusher	8.5 Fr	4.2 mm
SSP-10-X Length	Stent Pusher	10 Fr	4.2 mm
SSP-11.5-X Length	Stent Pusher	11.5 Fr	4.2 mm