

Biliary Dilator

Used to dilate papilla or biliary strictures through the scope



Technical Information

Tapering length - 3 cm.

Size - 5fr. - 6fr. - 7 fr. - 8.5 fr. - 9fr. - 10fr. - 11.5fr.

Radioopaque PTFE material .

Availability - Sterile - Single use only.

Individually packed.



Ordering Information - Mention the product code x quantity required.

Reference	Product Name	Tip Ø	Required Channel Ø	Working Length
S - BD- 5	Biliary Dilator 5 Fr.	2.5 Fr	3.2 mm	200 cm
S - BD- 6	Biliary Dilator 6 Fr.	3.5 Fr	3.2 mm	200 cm
S - BD- 7	Biliary Dilator 7 Fr.	4 Fr	3.2 mm	200 cm
S - BD- 8.5	Biliary Dilator 8.5 Fr.	5.5 Fr	3.2 mm	200 cm
S - BD- 9	Biliary Dilator 9 Fr.	6 Fr	3.2 mm	200 cm
S - BD- 10	Biliary Dilator 10 Fr.	6 Fr	4.2 mm	200 cm
S - BD- 11.5	Biliary Dilator 11.5 Fr.	8.5 Fr	4.2 mm	200 cm
S - BD- 5 - DBE	Biliary Dilator 5 Fr. For Double Balloon Enteroscope	2.5 Fr	3.2 mm	280 cm
S - BD- 6 - DBE	Biliary Dilator 6 Fr. For Double Balloon Enteroscope	3.5 Fr	3.2 mm	280 cm
S - BD- 7 - DBE	Biliary Dilator 7 Fr. For Double Balloon Enteroscope	4 Fr	3.2 mm	280 cm
S - BD- 8.5 - DBE	Biliary Dilator 8.5 Fr. For Double Balloon Enteroscope	5.5 Fr	3.2 mm	280 cm
S - BD- 9 - DBE	Biliary Dilator 9 Fr. For Double Balloon Enteroscope	6 Fr	3.2 mm	280 cm
S - BD- 10 - DBE	Biliary Dilator 10 Fr. For Double Balloon Enteroscope	6 Fr	4.2 mm	280 cm
S - BD- 11.5 - DBE	Biliary Dilator 11.5 Fr. For Double Balloon Enteroscope	8.5 Fr	4.2 mm	280 cm