

OmniView

CE₂₄₆₀

MULTIPLE BAND LIGATOR CARTRIDGE



For endoscopic ligation of esophageal varices, at or above the gastroesophageal junction.

This product is supplied non-sterile and intended for single use only.

To be used by or on order of a physician only.

Medelec
SYSTEMS PVT LTD

CONSISTING OF

The product consists of:

- Ligator cartridge with ultra-clear, wide view polycarbonate barrel preloaded with 6 high quality elastomeric bands.
CAUTION: The elastomeric bands are made of natural rubber latex.
- Loading catheter
- Irrigation needle

ORDERING INFORMATION

Minimum accessory channel 2.8 mm.

The product is available in various sizes to fit endoscopes ranging from 8.6 mm to 14 mm diameter.

<u>Ref</u>	<u>Description</u>
MBLC-6	1 cartridge with 6 bands, for endoscope diameter 9.5 ~ 13.0 mm
MBLC-6-XL	1 cartridge with 6 bands, for endoscope diameter 11.0 ~ 14.0 mm
MBLC-6-XS	1 cartridge with 6 bands, for endoscope diameter 8.6 ~ 9.2 mm

PACKAGING INFORMATION

One set per Tray.

Marketed By

Globax Marketing & Consultants

303, Third Floor, Somdutt Chamber - II, 9, Bhikaji Cama Place, New Delhi - 110066,

T: 011-49051596/26176463, M: 9811601439

Email: info@globax.in, Website: www.globax.in



Medelec
SYSTEMS PVT LTD