

ERCP Cannula



Used endoscopically to opacify the biliary duct with radioopaque contrast material for fluoroscopic view

Technical Information

Type - Tapered tip / Ball tip / Double lumen


Catheter tip Diameter - 4.5 Fr. / 5.5 Fr.

Recommended wire guide Diameter 0.035 inch.

Availability - Sterile For single use only.

Length- 200 cm & 280 cm

Individually packed.



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Tip Ø	Working Length	Required Channel Ø
S- ERC-T-SL	Single Lumen ERCP Cannula Tapered Tip	4.5 Fr	200 cm	2.8 mm
S- ERC-B-SL	Single Lumen ERCP Cannula Ball Tip	5.5 Fr	200 cm	2.8 mm
S- ERC-T-DL	Wire Guided ERCP Cannula Tapered Tip	5.5 Fr	200 cm	2.8 mm
S- ERC-T-SL-DBE	Single Lumen ERCP Cannula Tapered Tip For Double Balloon Enteroscope	4.5 Fr	280 cm	2.8 mm
S- ERC-B-SL-DBE	Single Lumen ERCP Cannula Ball Tip For Double Balloon Enteroscope	5.5 Fr	280 cm	2.8 mm
S- ERC-T-DL-DBE	Wire Guided ERCP Cannula Tapered Tip For Double Balloon Enteroscope	5.5 Fr	280 cm	2.8 mm