

Biliary Metal Stent



Used for stenosis or obstructed biliary tract caused
By benign or malignant changes

Technical Information

With Introducer System.

Radiotherapy compatible.

Expanded Stent Diameter. - 10mm.


Stent Length - 60mm / 80mm / 100mm

Uncovered stent.

Compatible with 0.038 inch wire guide.

Compatible with instrument channel size - 3.2mm

Availability - Sterile - For single use only.



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Stent Length	Stent Ø
S - MBS - 6	Biliary Metal Stent	60mm	10mm
S - MBS - 8	Biliary Metal Stent	80mm	10mm
S - MBS - 10	Biliary Metal Stent	100mm	10mm