

# Clipping Device

Used to mechanically bind blood vessels




## Technical Information

Working length 165/230cm

Sheath diameter 2.6 mm

Availability - sterile for single use only.

Clip opening diameter 15 mm



Ordering Information : Mention the product code x quantity required.

| Reference | Product Name    | Length | Sheath Diameter |
|-----------|-----------------|--------|-----------------|
| S-CD-16   | Clipping Device | 165cm  | 2.6mm           |
| S-CD-21   | Clipping Device | 230cm  | 2.6 mm          |