

Cystostome



To electro surgically cannulate the transgastric or transduodenal wall and into a pancreatic pseudocyst that is visibly bulging into the gastrointestinal tract

Technical Information

Size - 6 - 7 - 8 - 10 Fr.

Length - 180 cm

Availability - Sterile for single use only.



Ordering Information- - Mention product code x quantity required

Reference	Product Name	Size	Required Channel Ø	Wire Guide Compatibility
S - CYS - 6	Cystostome	6 Fr.	2.8mm	0.025"
S - CYS - 7	Cystostome	7 Fr.	2.8mm	0.025"
S - CYS - 8	Cystostome	8 Fr.	3.8mm	0.025"
S - CYS - 10	Cystostome	10 Fr.	3.8mm	0.035" for Single stenting
S - CYS - 10	Cystostome	10 Fr.	3.8mm	0.025" for Double stenting