

Esophageal Dilator

Used to dilate esophageal strictures over previously placed guide wire

Available Sizes - 5mm to 20mm Outer Diameter.
Other sizes can be produced on specific request.



Soft, Flexible and Radioopaque Tip & Radioopaque Junction.
Gradual Tapering.

Low friction lumen surface.

Compatible with 0.038 inch wire guide.

Reference No.:	Product Name
SED-SET-6-70	Esophageal Dilator set
SED-SET-6-100	Esophageal Dilator set
SED-SET-12-70	Esophageal Dilator set
SED-SET-12-100	Esophageal Dilator set
SED-SET-16-70	Esophageal Dilator set
SED-SET-16-100	Esophageal Dilator set

Supplied Non Sterile only

Set of 6 Pcs... Any 6 Pc from 5MM to 20MM

Set of 12 Pcs... Any 12 Pc from 5MM to 20MM

Set of 16 Pcs from 5MM to 20MM

Set includes 1pc Flexible spring tip Guide Wire 250 cm long, 1pc Cleaning Brush. Carrying case is available with order of set only.



Ordering Information- Mention the product code x quantity required.

Reference	Outer Ø	Length
SED-500	5 mm	100 cm
SED-600	6 mm	100 cm
SED-700	7 mm	100 cm
SED-800	8 mm	100 cm
SED-900	9 mm	100 cm
SED-1000	10 mm	100 cm
SED-1100	11 mm	100 cm
SED-1200	12 mm	100 cm
SED-1300	13 mm	100 cm
SED-1400	14 mm	100 cm
SED-1500	15 mm	100 cm
SED-1600	16 mm	100 cm
SED-1700	17 mm	100 cm
SED-1800	18 mm	100 cm
SED-1900	19 mm	100 cm
SED-2000	20 mm	100 cm

Reference	Outer Ø	Length
SED-570	5 mm	70 cm
SED-670	6 mm	70 cm
SED-770	7 mm	70 cm
SED-870	8 mm	70 cm
SED-970	9 mm	70 cm
SED-1070	10 mm	70 cm
SED-1170	11 mm	70 cm
SED-1270	12 mm	70 cm
SED-1370	13 mm	70 cm
SED-1470	14 mm	70 cm
SED-1570	15 mm	70 cm
SED-1670	16 mm	70 cm
SED-1770	17 mm	70 cm
SED-1870	18 mm	70 cm
SED-1970	19 mm	70 cm
SED-2070	20 mm	70 cm