

# Esophageal Dilator



Used to dilate esophageal strictures over previously placed guide wire

Available Sizes - 5mm to 20mm Outer Diameter.

Other sizes can be produced on specific request.

Soft, Flexible and Radioopaque Tip & Radioopaque Junction.

Gradual Tapering.

Low friction lumen surface.

Compatible with 0.038 inch wire guide.



Dilator Size - Diameter - 5 mm. to 20 mm. Length 100 cm. , 70 cm.

Availability - Sterile for single use only

**Ordering Information- Mention the product code x quantity required.**

| Reference  | Outer Ø | Length |
|------------|---------|--------|
| SED-500-S  | 5 mm    | 100 cm |
| SED-600-S  | 6 mm    | 100 cm |
| SED-700-S  | 7 mm    | 100 cm |
| SED-800-S  | 8 mm    | 100 cm |
| SED-900-S  | 9 mm    | 100 cm |
| SED-1000-S | 10 mm   | 100 cm |
| SED-1100-S | 11 mm   | 100 cm |
| SED-1200-S | 12 mm   | 100 cm |
| SED-1300-S | 13 mm   | 100 cm |
| SED-1400-S | 14 mm   | 100 cm |
| SED-1500-S | 15 mm   | 100 cm |
| SED-1600-S | 16 mm   | 100 cm |
| SED-1700-S | 17 mm   | 100 cm |
| SED-1800-S | 18 mm   | 100 cm |
| SED-1900-S | 19 mm   | 100 cm |
| SED-2000-S | 20 mm   | 100 cm |

| Reference  | Outer Ø | Length |
|------------|---------|--------|
| SED-570-S  | 5 mm    | 70 cm  |
| SED-670-S  | 6 mm    | 70 cm  |
| SED-770-S  | 7 mm    | 70 cm  |
| SED-870-S  | 8 mm    | 70 cm  |
| SED-970-S  | 9 mm    | 70 cm  |
| SED-1070-S | 10 mm   | 70 cm  |
| SED-1170-S | 11 mm   | 70 cm  |
| SED-1270-S | 12 mm   | 70 cm  |
| SED-1370-S | 13 mm   | 70 cm  |
| SED-1470-S | 14 mm   | 70 cm  |
| SED-1570-S | 15 mm   | 70 cm  |
| SED-1670-S | 16 mm   | 70 cm  |
| SED-1770-S | 17 mm   | 70 cm  |
| SED-1870-S | 18 mm   | 70 cm  |
| SED-1970-S | 19 mm   | 70 cm  |
| SED-2070-S | 20 mm   | 70 cm  |