

OmniTip ^{UltraHD}

CE₂₄₆₀

Echogenic Ultra Sound Needle




An Endoscopic Ultra Sound Needle system used to sample targeted submucosal gastrointestinal lesions through the accessory channel of an ultrasound endoscope.

This device is supplied sterile and is disposable – intended for single use only.

To be used by, or on order of, a physician only.

Medelec
SYSTEMS PVT LTD

- Needle length extension adjustable up to 8 cm using the safety lock slider.
- Sheath length extendable up to 5 cm to fit a full range of endoscopes.
- Enhanced echogenic needle tip for clear ultrasonic visualization.
- Ergonomically contoured handle for better control and stability.
- Luer lock for secure and stable attachment to the endoscope.



ORDERING INFORMATION

Minimum accessory channel 2.0 mm. 10 cc vacuum locking syringe included.

Ref No.	Needle Description	Sheath Size
OTU 19-A	19 gauge needle with ball tip stylet	5.2 FR
OTU 22-A	22 gauge needle with ball tip stylet	5.2 FR
OTU 24-A	24 gauge needle with ball tip stylet	5.2 FR
OTU 25-A	25 gauge needle with ball tip stylet	5.2 FR

PACKAGING INFORMATION

One set of EUS needle and Vacuum syringe per Tray. Two trays per Box.
Five Boxes per Carton.

Marketed By

Globax Marketing & Consultants
303, Third Floor, Somdutt Chamber- II,9, Bhikaji Cama Place, New Delhi - 110066,
Tel: 011-049051596/26176463, M: 9811601439
Email: info@globax.in, Website: www.globax.in



Medelec
SYSTEMS PVT LTD

EUS NEEDLE



Used for aspiration biopsy to diagnose and stage GI lesions, the tip provides sampling by employing ultrasound technology to guide the needle

HANDLE IS EASY TO DETACH TO REUSE AND IS ECONOMICAL, SAFE & CONTAMINATION FREE

Technical Information



Working Length - 140 cm

Needle Tip - 19 & 22 swg

Handle is Autoclavable up to 10 times at 140°C if it is in working condition.

Sheath/Needle & Stylet is available in Sterile pack

Type - FNA, FNB



Ordering Information- - Mention product code x quantity required

Reference	Product Name	Type	Sheath Ø	Needle Gauge	Length
S-EUS-S-19-A	EUS Needle	FNA	2.30 MM	19 SWG	140 cm
S-EUS-S-19-B	EUS Needle	FNB	2.30 MM	19SWG	140 cm
S-EUS-S-22-A	EUS Needle	FNA	2.00 MM	22 SWG	140 cm
S-EUS-S-22-B	EUS Needle	FNB	2.00 MM	22 SWG	140 cm
S-EUS-LA-19-A	Sheath with Needle & Stylet	FNA	2.30 MM	19 SWG	140 cm
S-EUS-LA-19-B	Sheath with Needle & Stylet	FNB	2.30 MM	19 SWG	140 cm
S-EUS-LA-22-A	Sheath with Needle & Stylet	FNA	2.00 MM	22 SWG	140 cm
S-EUS-LA-22-B	Sheath with Needle & Stylet	FNB	2.00 MM	22 SWG	140 cm