

Cytology Brush



Used to collect cells from the upper and lower gastrointestinal tracts

Technical Information

Length-170cm / 230cm & 325cm

Wire Guide Compatibility-0.035"

Availability- Individually packed. Sterile For single use only.



Ordering Information : Mention the product code x quantity required.

Reference	Product Name	Required Channel Ø	Length
S - CY-170	Cytology Brush	2.8 mm	170 cm
S - CY-230	Cytology Brush	2.8 mm	230 cm
S - CY-DL	Wire Guided Cytology Brush	2.8 mm	230 cm
S - CY - S	Cytology Brush For Spy	1.2 mm	325 cm