


Balloon Inflation Device

Technical Information :

Reusable, Non-Sterile

For Achalasia Balloon Dilator



PNEUMATIC

Technical Information :

Controlled push of air/water.

0 to 30 atm indicator guage meter

20ml syringe with 3 way stop-cock

Non-Sterile for single use only

For Esophageal, Biliary Balloon Dilator &
3 Stage Balloon Dilator



MICRO

Ordering Information- Mention the product code quantity required.

Reference	Product Name
SBID-IN	Balloon Inflation Device - Pneumatic
SBID	Balloon Inflation Device - Micro