

# Sclerotherapy Needle



**Injector needles- sclero therapy needle**  
**used for endoscopic injection into esophageal varices above O.G. junction**


Used for Endoscopic injection to introduce a sclerosing agent of vasoconstrictor into selected sites to control actual or potential bleeding lesions. The injection of saline to aid in Endoscopic Mucosal Resection (EMR), Polypectomy procedures and to control non-variceal haemorrhage.

## Technical specifications -

**Needle** - Short and Sharp Bevel, Available in 19/20/21/22/23/24/25 SWG.

**Working length** - 170 cm. / 230 cm. / 280 cm

**Availability** -Sterile For Single use only.  
 Specially Developed Needle To Provide Secure Firm Grip With Inner Tube & Prevent Possible Leakage From The Joint Of Inner Tube & Needle.  
 Specially Developed Needle Give "JET" Pressure To Inject The Drug for 23/24 SWG.



**Ordering Information- Mention the product code x quantity required.**

| Reference    | Product Name                   | Needle Gauge | Working Length | Required Channel Ø |
|--------------|--------------------------------|--------------|----------------|--------------------|
| SGN-2319     | Sclerotherapy Needle           | 19G          | 170 cm         | 2.8 mm.            |
| SGN-2319-C   | Sclerotherapy Needle           | 19G          | 230 cm         | 2.8 mm.            |
| SGN-2320     | Sclerotherapy Needle           | 20G          | 170 cm         | 2.8 mm.            |
| SGN-2320-C   | Sclerotherapy Needle           | 20G          | 230 cm         | 2.8 mm.            |
| SGN-2321     | Sclerotherapy Needle           | 21G          | 170 cm         | 2.8 mm.            |
| SGN-2321-C   | Sclerotherapy Needle           | 21G          | 230 cm         | 2.8 mm.            |
| SGN-2321-M   | Sclerotherapy Needle Metal Tip | 21G          | 170 cm         | 2.8 mm.            |
| SGN-2321-C-M | Sclerotherapy Needle Metal Tip | 21G          | 230 cm         | 2.8 mm.            |
| SSN-2322     | Sclerotherapy Needle           | 22G          | 170 cm         | 2.8 mm.            |
| SSN-2322-C   | Sclerotherapy Needle           | 22G          | 230 cm         | 2.8 mm.            |
| SSN-2323     | Sclerotherapy Needle           | 23G          | 170 cm         | 2.8 mm.            |
| SSN-2323-C   | Sclerotherapy Needle           | 23G          | 230 cm         | 2.8 mm.            |
| SSN-2323-M   | Sclerotherapy Needle Metal Tip | 23G          | 170 cm         | 2.8 mm.            |
| SSN-2323-C-M | Sclerotherapy Needle Metal Tip | 23G          | 230 cm         | 2.8 mm.            |
| SSN-2324     | Sclerotherapy Needle           | 24G          | 170 cm         | 2.8 mm.            |
| SSN-2324-C   | Sclerotherapy Needle           | 24G          | 230 cm         | 2.8 mm.            |
| SSN-2325     | Sclerotherapy Needle           | 25G          | 170 cm         | 2.8 mm.            |
| SSN-2325-C   | Sclerotherapy Needle           | 25G          | 230 cm         | 2.8 mm.            |
| SSN-2324-DBE | Sclerotherapy Needle           | 24G          | 280 cm         | 2.8 mm.            |