

Foreign Body Holding Forcep



Used to retrieve foreign body from G.I. tract.

Shark Tooth Forcep.

Alligator Forcep.

Shark Tooth Alligator Forcep.

Features :


Can be used with ERCP Scope to Retrieve stent.(Shark Tooth Forcep Only).

Supplied Without Sterile and it is Reusable.

Specially produced outer sheath for easy elevation through elevator. (Shark Tooth Forcep Only)

Firm grip & easy movement.

Rotatable graspers are available



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Required Channel Ø	Jaws Specification	Working Length
SSTF-230	Shark Tooth Forcep	2.8 mm	Rat Tooth Jaws	230 cm
SSTF-170	Shark Tooth Forcep	2.8 mm	Rat Tooth Jaws	170 cm
SATF-230	Alligator Forcep	2.8 mm	Alligator Jaws	230 cm
SATF-170	Alligator Forcep	2.8 mm	Alligator Jaws	170 cm
SSTF-230-R	Shark Tooth Forcep Rotatable	2.8 mm	Rat Tooth Jaws	230 cm
SSTF-170-R	Shark Tooth Forcep Rotatable	2.8 mm	Rat Tooth Jaws	170 cm
SATF-230-R	Alligator Forcep Rotatable	2.8 mm	Alligator Jaws	230 cm
SATF-170-R	Alligator Forcep Rotatable	2.8 mm	Alligator Jaws	170 cm
SSAF-230	Shark Tooth Alligator Combo	2.8 mm	Combo	230 cm
SSAF-170	Shark Tooth Alligator Combo	2.8 mm	Combo	170 cm
SSAF-230-R	Shark Tooth Alligator Combo Rotatable	2.8 mm	Combo	230 cm
SSAF-170-R	Shark Tooth Alligator Combo Rotatable	2.8 mm	Combo	170 cm

“Shaili “ Range of Foreign Body Holder : 2 prong , 3 Prong , Retrieval Basket, Snare, Shark Tooth Forcep, Alligator Forcep, Monofilament Memory Basket, Loop Basket.