



8 Endoscope Storage & Drying Cabinet

- Stainless Steel 304 Grade Body
- Total Microbiological Safety
- High Pressure Blow Fans
- Rotatable Hanger for easy use
- Air flow through HEPA Filter
- LED Lights with Controller
- Toughened Glass Door with lock

Manufactured By:

Endotech



Technical Specifications

Endoslide 8

8 Endoscope Storage Cabinets

Technical Details for Models (Body Stainless steel grade ss304):

1. Fresh Air Fan with HEPA Filter
2. Exhaust Fan
3. Circulation Fan for Drying
4. LED Lights - Door Open & Close
5. Timer for Fans & LED
6. Door - toughened glass Locking with key
7. for 8 scope vertical rotatable hanger